


XpressRNA Viral Kit





- 100% Sensitivity and Specificity
- More than 95% effective
- Compatible for Manual and Automation
- Suitable for ZYBIO, KingFisher, Himedia & ALTA
- For HCV, HPV, HBV, HIV, H1N1, HLA typing & TB

- Sample type : Saliva, Sputum, Swabs, Serum, Plasma, Whole Blood and liquid
- Sample volume : 200µL (VTM/MTM/UTM/Serum/Plasma/Blood/HPV Media)
- Elution volume : 50µL
- Throughput : 32 - 96 samples per run
- Extraction time : 15-30 Mins*

Patented Magnetic Nanoparticles

- Iron oxide MNP: Biocompatible and FDA approved
- Average size 10-12 nm
- Superparamagnetic - Particles gets magnetized on application of magnetic field, while on removal it behaves like non-magnetic particles
- High surface to volume ratio



ORDERING INFORMATION

CAT No	DESCRIPTION	PACK SIZE
MG20VRNA50	XpressRNA Viral Kit	50 Tests
MG20VRNA250	XpressRNA Viral Kit	250 Tests